
# Molift EvoSling Comfort MediumBack

**molift**<sup>®</sup>  
by Etac


User Manual - Bruksanvisning - Brukermanual - Brugervejledning - Käyttö-  
pas - Bedienungsanleitung - Gebruikershandleiding - Manuel d'utilisation  
- Manuale per l'utente - Manual de usuario



BM28699 Rev. B 2021-02-22

MD

CE

## Contents

User manual - English - 4
Bruksanvisning - Svenska - 8
Brukermanual - Norsk - 12
Brugsvejledning - Dansk - 16
Käyttöohje - Suomi - 20
Gebrauchsanweisung - Deutsch - 24
Handleiding - Nederlands - 28
Manual d'utilisation - Français - 32
Manual utente - Italiano - 36
Manual de usuario - Español - 40

## General

### Declaration of conformity

The Molift EvoSling Comfort MediumBack and related accessories described in this operator manual are CE marked in accordance with EU Council Directive MDR (EU) 2017/745 concerning medical devices, class 1, and has been tested and approved by a third party according to standard EN ISO 10535:2006.

**Conditions for Use:** Lift and transfer of a person will always pose a certain risk and only trained personnel should use the equipment and accessories covered by this user manual.

**Warranty:** 2-year warranty against defects in workmanship and materials of our products. For Terms and conditions,

see [www.etac.com](http://www.etac.com)

Visit [www.etac.com](http://www.etac.com) to download the current documentation.

## Important

This User Manual contains important safety instructions and information regarding the use of the sling and accessories. In this manual the user is the person being hoisted. The assistant is the person operating the hoist.

This symbol indicates important information related to safety. Follow these instructions carefully.

**Read User Manual before use!**  
It is important to fully understand the content of the user manual before attempting to use the equipment. Read the manual for both hoist and sling

**General Safety Precautions:** Molift EvoSling MediumBack is designed for use with a 2-point suspension but also works with a 4-point suspension. Only use accessories and slings that are adjusted to fit the user, type of disability, size, weight and type of transfer. It is important that the sling has been tested with the individual user and for the intended hoisting situation. Make a decision on whether one or more assistants are required. Plan the hoisting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.

Assess the risk and take notes.

You as the assistant are responsible for the safety of the user!

Molift slings should only be used to hoist people. Never use the sling to lift or move objects of any kind.


If maximum load (SWL) differs between hoist, suspension and body support unit, then the lowest maximum load shall always be used

The sling should not be stored in direct sunlight.

## Labels and symbols

1. Main label
2. Size label
3. Batch and serial number label w. GS1-128 barcode.
4. Product label
5. Periodic inspection label/ Name tag

### Symbols:

Manufacturer

Date of manufacture.  
YYYY-MM-DD  
(year/month/date)

CE marked

Refer to user manual

This way up, this side is outside

Max user weight

Medical Device

### Wash symbols:

Slings can be washed at temperatures ranging from 60 up to 85° C

Tumble drying, max 60° C.

Do not iron

Do not dry clean


Do not bleach


## About EvoSling Comfort MediumBack

Etac supplies a wide selection of slings for different types of transfers. Molift EvoSling Comfort MediumBack is a sling designed to give a comfortable transfer and to remain under the user after transfer into a sitting position, ie, in a wheelchair. The sling does not have split leg support and has a shoulder height back which gives support for the seat area and entire back. The sling can be used to hoist from sitting position in bed to sitting position out of bed. Approved for users with weight up to 300kg.

Molift EvoSling Comfort MediumBack Net is available in sizes XS – XXL, made of net polyester. The sling has a centre mark line to facilitate application of the sling. The sling has a guide handle for positioning of the user when hoisting. The user must have their arms inside the sling. The sling can only be applied when the user is in a lying position. Molift EvoSling is developed to be used with a 2-point suspension but the sling also works with a 4-point suspension. See the combination list for the correct sling and suspension.

**Technical data**

**Safe Working Load (SWL):**  
XXS-XL: 300 kg (660 lbs)

**Material:**  
Polyester, Polyethylene

Size	A	B
XS	77	83
S	90	86
M	102	94
L	115	113
XL	118	120
XXL	120	126


**Expected Lifetime:**

The product has an expected service life of 1 to 5 years under normal use. The service life of the product varies depending on usage frequency, materials, loads and how often it is washed. Washing at higher temperatures wears out the material faster.

**Lifting Strap Ladder loops**

The ladder loops on the lifting band give different options for mounting the sling to the suspension. Use the ladder that makes the sling most comfortable for the user.

The ladder loops on the straps have different colours so the strap can be easily fastened with equal length on left and right side. Distance between each ladder loop (A) is 9cm.

**Combination list**

 2 Point Suspension	Sling size		
	XS-S	M-L	XL-XXL
Small 340mm	✓	✓	
Medium 440mm	✓	✓	✓
Large 540mm		✓	✓

 4 Point Suspension	Sling size		
	XS-S	M-L	XL-XXL
Small 340mm	✓	✓	
Medium 440mm	✓	✓	✓
Large 540mm		✓	✓

**Before use / Daily check**

 Before using the sling for the first time it should be marked with the date of first time use on the periodic inspection label.

Inspection to be performed daily or before use:

- Make sure the user has the required ability for the sling
- Make sure the sling is suitable for the hoist/ suspension to be used.
- Make sure sling does not have visible damage or frays.


**Periodic inspection**

Periodic inspection of the non-rigid body support unit should be undertaken at the time intervals stated by the manufacturer, but at least every 6 months. More frequent inspections may be required where a non-rigid body support unit is used or cleaned more frequently than normal

 Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings

 Make sure not to damage or remove labels when cleaning.

**Size guide**


Measurements in the table should be used as a reference. Correct size depends on the user's weight, function and body shape. Measure the waist size and back height on user in a sitting position. Try the sling with the most correct measurements to make sure it fits.

Size	User weight - kg (lbs)	Width (W) - cm	Height (H) - cm
XS	17-25 (37-55)	31-38	60-70
S	25-50 (55-110)	35-42	70-80
M	45-95 (99-210)	37-47	80-95
L	90-165 (198-353)	40-51	95-100
XL	160-240 (353-529)	45-70	100-105
XXL	230-300 (507-660)	50-75	100-120


**Sling checkpoints**

1. Ensure the suspension is away from users head
2. The leg support is under the users thigh without creases

## Transfer

When moving the user, stand to the side of the person you are lifting. Make sure that arms and legs do not obstruct the seat, bed, etc.

Never lift the user higher than necessary to carry out a lift. Remember that wheels on a mobile lift must NOT be locked.


1. Check that the sling is correctly fitted around the user and that the strap loops are correctly fitted to the suspension hooks.
2. Start lifting until the lifting straps are stretched without lifting the user. Ensure that all four loops of the sling are securely fastened to avoid the user slipping or falling out of the sling.
3. Hoist the user, and perform transfer.

Be careful during movement, the user may swing during turns, stops and starts. Be careful when manoeuvring close to furniture and similar to ensure the suspended user does not hit anything.

Avoid deep pile carpets, high thresholds, uneven surfaces or other obstacles that may block the castors. The lifter may become unstable if forced over such obstacles increasing the risk of tipping over.

User must be facing the hoist when lifting or lowering from or into a chair.


Never leave a user unattended while being hoisted.


The hoist should not be used to lift or move users on sloping surfaces

## Lowering into sitting


When transferring the user into a sitting position, the following techniques can be used for lowering the user as far back into the seat as possible (Use A, B or C):

A Push gently on the users knees to steer the user in towards the backrest of the chair.


B Tilt the seat backwards.

C Position the user with the help of the steering handle

Do not lift with the guide handles! Excessive force on the handles can cause the sling to tear.

The sling is designed to remain under the user after transfer into a sitting position.


If the sling has been displaced under the user there may be a risk of them falling out of the sling when being lifted again.

Make sure the sling is correctly positioned on a user before lifting. Transfer the user into lying position to remove the sling.


## Lying position - Sling application


It is only recommended to use MediumBack slings from a lying position if the user can be transferred in a seated position.

1. Turn the user towards you so that the user is lying in a safe position on their side. Place the sling centered on the users back with the bottom edge placed above the knee under the users thigh. Fold the sling so it is possible to pull out the folded side after turning the user over on the other side.

2. Pull out the folded sling. Make sure the sling is placed in the centre of the users back.

3. Make sure the lifting straps are the same length on both sides and that the sling is placed under the users thighs without any folds. Raise the user into sitting position before lifting.

Make sure the sling is placed well under the thighs close to the knee. User can slip or fall out of sling if it is not placed under the entire length.


Connect all 4 lifting straps to suspension. Make sure the lifting straps are fastened with equal length on right and left side. Upper straps and leg straps can be mounted with different colours on ladder loop to make sure sling is comfortable and safe.

Do not cross the leg support. Leg straps will squeeze the user and might cause pain.

## Lying position - Removing the sling

Lower the user into lying position. Release the lifting straps from the suspension.

1. Turn the patient towards you so that the user is lying in a safe position on their side. Fold the sling in the middle and position it behind the user's back.

2. Turn the user over on the opposite side and pull the sling away. Support the user with one hand if necessary.

Remove the sling gently to ensure the user does not roll over and fall off the bed.

## Recycling

Recycling instructions is available at [Etac.com](http://Etac.com)

## Troubleshooting

Symptom	Possible Cause/Action
The sling does not fit the user.	Wrong sling size. Try another size.
	The sling is too small when the suspension is too close to the user's face. Try another size.
All the straps do not lift at the same time.	If the sling straps lift unevenly, it means that the sling can slide up at the back. Make sure that the sling is positioned correctly under the thighs.
User is sliding through/out of the sling.	The sling is not connected correctly according to the users prescription.
	The user is wearing "slippery" clothes. Change or remove some of the clothing
	Sling has been displaced under user when sitting on sling. Adjust sling, and makes sure it is correctly positioned on user.

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Etac A/S