Molift UnoSling HighBack

molift[®]
by Etac

User manual

BM29599 Rev D 2021-05-10

MD

CE

Index

Labels and symbols.....	3
User manual – English.....	4
Bruksanvisning – Svenska.....	6
Brukermanual – Norsk.....	8
Brugervejledning – Dansk.....	10
Käyttöohje – Suomi.....	12
Gebrauchsanweisung – Deutsch.....	14
Handleiding – Nederlands.....	16
Manual d'utilisation – Français.....	18
Manual utente – Italiano.....	20
Manual de usuario – Español.....	22

Labels and symbols/ Etiketter och symboler/ Etiketter og symboler/ Mærkater og symboler/ Tarrat ja symbolit/ Beschriftungen und Symbole/ Labels en symbolen/ Étiquettes et symboles/ Etichette e simboli/ Símbolos y etiquetas:

Molift UnoSling HighBack
A patient specific, disposable sling
Art. no: xxxxxxxx
DK HMI: xxxxxx

Size: XX

Material: Polypropylene, Polyester

Barcode according to GS1-128
(01)xxxxxxxxxxxx(11)xxxx(21)xxxxxx
(01)EAN number (11)YYMMDD(21)Serial number

UP/UPP
OUT/UT

YYYY-MM-DD

300 kg 66 lbs

Date of first use

Year	2021	2022	2023	2024	2025	2026	2028	2028
Month	01	02	03	04	05	06	07	08
	09	10	11	12				

Etac A/S
Parallevej 3
DK-8751 Gedved

MD

CE

DO NOT WASH DO NOT BLEACH DO NOT TUMBLE DRY DO NOT IRON DO NOT DRY CLEAN

Name: _____

Symbols:

 EN - Manufacturer
SV - Tilverkare
NO - Produsent
DA - Producent
FI - Valmistaja
DE - Hersteller
NL - Fabrikant
FR - Fabricant
IT - Produttore
ES - Fabricante

 EN- This Way up, this side is outside
SV - Denna sida upp, denna sida utåt
NO - Denne siden opp, denne siden ut
DA - Denne side op, denne side er ydersiden
FI - Tämä puoli ylöspäin, tämä puoli ulospäin
DE - Hier oben, diese Seite ist außen
NL - Deze kant boven, dit is de buitenkant
FR - Vers le haut, côté extérieur
IT - Lato alto, lato esterno
ES - Este lado hacia arriba, este lado fuera

Wash symbols:

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
Molift UnoSling HighBack

The Molift UnoSling is a consumable for single patient use.

The sling has a high back with head support and split leg support with a medium opening. It provides a reclined sitting position with support for the whole body. The loops on the lifting straps give two different height settings. Use the loop that makes the sling most comfortable for the user.

The UnoSling HighBack is recommended for lifting from a sitting or lying position to a sitting or lying position. The sling is not suitable for use in conjunction with bathing or showering, or for users who have had a leg amputated.

The Molift UnoSling HighBack sling is a consumable for single patient use, it has a limited service life.


Find the correct size and write the user's name and the date on the product label before using the product.

The UnoSling HighBack is manufactured from non-woven material. Ensure that no pointed or sharp objects come into contact with the sling. Avoid pulling hard on the sling when positioning the user in a bed or wheelchair.

The UnoSling HighBack should be taken out of use:

- when it has become dirty or damp, or if there is a risk that the sling has been exposed to contamination;
- if it has been damaged;
- when the user no longer requires it.

The sling should be disposed of in accordance with applicable procedures for the health care facility.


Things to remember

- Read through the manual before starting to use the hoist.
- It is important to understand the contents of the user manual before using the equipment. Read the manual for both hoist and sling.
- Make a decision on whether one or more assistants are required. Plan the lifting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.
- Lifting and transferring a person always involves a certain risk. Only informed personnel should use the equipment described in this user manual. If the maximum load (SWL) differs between the suspension and the sling, the lowest maximum load should always be used.
- Carry out a risk assessment and document it. As the assistant, you are responsible for the safety of the user.
- Never leave a user unattended in a lifting situation.
- Make sure that the wheels on the wheelchair, bed, etc. are locked during the movement.
- Remove the sling carefully to ensure that the user does not fall forwards.


The following inspection should be carried out daily or before use

- The user must have the bodily functions required for the use of the sling.
- Make sure the sling is suitable for the hoist/suspension to be used.
- Check that the sling has no visible damage or faults.

Sling checkpoints

Checkpoint	Troubleshooting
1 Safety clearance to hoist – minimises the risk that the user will hit their head.	If the sling is too small, the head may come too close to the suspension.
2 The lower edge of the sling should be at the level of the coccyx, and should enclose the buttocks.	If the sling is too large, the opening at the buttocks will also be too large, and there may be a risk that the user will slide out. If the sling is too small, it may cut into the crotch.
3 Leg supports placed under the user's thighs without folds.	If the leg supports are not stretched tense, the sling may be pulled out of position before the lift, which can result in the risk of the user sliding out.
4 Connection of leg supports.	Crossed leg supports minimise the risk of sliding out or sliding forwards in the sling.
5 The lifting straps are stretched simultaneously when the lift is started. See 'Lifting checkpoints' below.	If the straps are stretched unequally, the sling may be pulled out of position, and there may be a risk of the user sliding out.


Lifting checkpoints

- 1** Check that the sling is correctly fitted around the user and that the lifting loops are correctly fitted to the suspension hooks.
- 2** Start lifting. Ensure that the lifting straps are stretched tense, and pause before lifting the user. Check that all lifting loops are attached correctly to the suspension before lifting the user.
- 3** Lift the user and carry out the planned transfer.

Lifting from a sitting position


- 1 Ask the user to lean forwards, and support them if necessary. Slide the sling down with the lower edge to the user's coccyx. The upper part of the sling supports the user's head.
- 2 Pull the leg supports forwards and position them under the user's thighs.
- 3 Make sure the lifting straps are the same length on both sides, and that the leg supports are positioned under the user's thighs without any folds. Cross the leg loops by pulling one of the loops through the other and connect all 4 lifting straps to the suspension. Make sure the user is in a comfortable and secure position in the sling. The user should be facing the hoist when being lifted from or to a sitting position in a chair or wheelchair.

Lifting to a sitting position


When the user is lowered into a sitting position, the following methods can be used so that the user is set down as far back in the seat as possible (A or B):

- A** Push gently on the user's knees to steer them in towards the backrest of the chair.
- B** Tilt the seat backwards slightly as the user is lowered into it. Remove the leg supports by folding them under the sling and the user's thighs, and pulling them out. Support the user if necessary, and then remove the sling by carefully pulling it upwards.

Lifting from a lying position


- 1 Turn the user towards you so that they are lying in a stable position on their side. Position the sling with the bottom edge in line with the user's coccyx. Make sure the sling is centred on the user's back, and that the top of the sling is supporting the user's head. Fold the sling so that you can pull out the folded side after turning the user over onto their back.

- 2 Pull the leg supports forwards and position them under the user's thighs. Cross the lifting straps of the leg supports by pulling one of the straps through the other, and connect all four lifting straps to the suspension. Raise the user closer to a sitting position before lifting if possible.

Lifting into a lying position

Lower the user into a sitting or lying position in the bed. Pull out the leg supports towards the sides. Turn the user towards you so that they are lying in a stable position on their side. Roll the sling behind the user's back, turn the user onto their back and carefully pull the sling away.

The Molift UnoSling HighBack and related accessories, as described in this manual, are CE labelled and comply with the Council of Europe's directive on medical devices MDR (EU) 2017/745, class 1, and have been tested and approved by a third party in accordance with the EN ISO 10535:2006 standard.

- This symbol indicates especially important safety information. Follow these instructions carefully.
- Do not wash. If the sling is washed, this symbol disappears and another symbol appears. The sling should then be disposed of.
- Do not use.

Safe Working Load (SWL): S-L: 300kg (660lbs) XL-XXL:350kg (771lbs)
Material: Polypropylene, polyester


2-point suspension 	Size of sling		
	S	M-L	XL-XXL
Small 340 mm	✓	✓	✗
Medium 440 mm	✓	✓	✓
Large 540 mm	✗	✓	✓

4-point suspension 	Size of sling		
	S	M-L	XL-XXL
Small 340 mm	✓	✓	✗
Medium 440 mm	✓	✓	✓
Large 540 mm	✗	✓	✓
X-Large 640 mm	✗	✓	✓

Item no.	Sling model	Size
3060020	Molift UnoSling HighBack Single-patient sling, consumer item	S
3060030	Molift UnoSling HighBack Single-patient sling, consumer item	M
3060040	Molift UnoSling HighBack Single-patient sling, consumer item	L
3060050	Molift UnoSling HighBack Single-patient sling, consumer item	XL
3060060	Molift UnoSling HighBack Single-patient sling, consumer item	XXL

Size	User weight, kg (lbs)	Width (W), cm	Height (H), cm
S	25–50 (55–110)	35–42	70–80
M	45–95 (99–210)	37–47	80–95
L	90–165 (198–353)	40–51	95–100
XL	160–350 (353–771)	45–70	100–105
XXL	230–350 (507–771)	50–75	100–120

Visit www.etac.com to download documentation to ensure that you have the latest version.

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