
Molift RgoSling MediumBack PLUS

molift[®]
a part of Etac

User manual

BM41399 Rev. A 2020-01-29

CE

Contents

User manual - English	- 4
Bruksanvisning - Svenska	- 8
Brukermanual - Norsk	- 12
Brugsvejledning - Dansk	- 16
Käyttöohje - Suomi	- 20
Gebrauchsanweisung - Deutsch	- 24
Handleiding - Nederlands	- 28
Manual d`utilisation - Francais	- 32
Manual utente - Italiano	- 36
Manual de usuario - Español	- 40

General

Declaration of conformity

The Molift RgoSling MediumBack Plus and related accessories described in this user manual are CE marked in accordance with EU Council Directive 93/42/EEC concerning medical devices, class 1, and has been tested and approved by a third party according to standard EN ISO 10535:2006.

Conditions for Use: Lift and transfer of a person will always pose a certain risk and only trained personnel should use the equipment and accessories covered by this user manual.

Warranty: 2-year warranty against defects in workmanship and materials of our products. For Terms and conditions, visit www.etac.com

Visit www.etac.com to download the current documentation.

Important

This User Manual contains important safety instructions and information regarding the use of the sling and accessories. In this manual the user is the person being hoisted. The assistant is the person operating the hoist.

This symbol indicates important information related to safety. Follow these instructions carefully.

Read User Manual before use!
It is important to fully understand the content of the user manual before attempting to use the equipment. Read the manual for both hoist and sling

General Safety Precautions: Molift RgoSling MediumBack Plus is designed only for use with a 4-point suspension. Only use accessories and slings that are adjusted to fit the user, type of disability, size, weight and type of transfer. It is important that the sling has been tested with the individual user and for the intended hoisting situation. Make a decision on whether one or more assistants are required. Plan the hoisting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.

Assess the risk and take notes.

You as the assistant are responsible for the safety of the user!

Molift slings should only be used to hoist people. Never use the sling to lift or move objects of any kind.


If maximum load (SWL) differs between hoist, suspension and body support unit, then the lowest maximum load shall always be used.

The sling should not be stored in direct sunlight.

Labels and symbols

1. Main label
2. Size label
3. Batch label
4. Product label with barcode that contains serial number
5. Product label with barcode that contains EAN number.
6. Periodic inspection label / Name tag

Symbols:

Date of manufacture. YYYY-MM (year/month)

CE marked

Refer to user manual

This way up, this side is outside

Max user weight

Wash symbols:

Slings can be washed at temperatures ranging from 60 up to 85° C.

Tumble drying, max 60° C.

Do not iron

Do not dry clean

Do not bleach


About RgoSling MediumBack Plus

Etac supplies a wide selection of slings for different types of transfers. Molift RgoSling MediumBack Plus has a shoulder-height back, split legsupport with two legstraps on each legsupport. The sling can be used for lifting to and from a sitting position.


Approved for users weighing up to 500kg.

Molift RgoSling MediumBack Plus is available in sizes L and XL, in durable polyester fabric with padded legsupport.

Molift RgoSling MediumBack Plus has been developed to be used with a 4-point steel sling bar. See the combination list for the correct sling and sling bar.

Technical data

Safe Working Load (SWL):

L-XL: 500 kg (1102 lbs)

Material:

Polyester, Polyethylene

Size	A	B
L	65	91
XL	68	107


Expected Lifetime:

The product has an expected service life of 1 to 5 years under normal use. The service life of the product varies depending on usage frequency, materials, loads and how often it is washed. Washing at higher temperatures wears out the material faster.


Lifting Strap Ladder loops

The ladder loops on the lifting strap give different options for mounting the sling to the suspension. Use the ladder that makes the sling most comfortable for the user.

The ladder loops on the straps have different colours so the strap can be easily fastened with equal length on left and right side. Distance between each ladder loop (A) is 9cm.

Combination list

 4-point steel sling bar	Sling size	
	L	XL
Art. no.: 2510435 SWL 500 kg	✓	✓

 4-point sling bar	Sling size	
	L	XL
Large 550mm Art. no.: 2530212 SWL 350 kg	✓	✓
X-Large 650mm Art. no.: 2530213 SWL 350 kg	✓	✓

Before use / Daily check

Before using the sling for the first time it should be marked with the date of the first time use on periodic inspection label.

Inspection to be performed daily or before use:

- Make sure the user has the required ability for the sling
- Make sure the sling is suitable for the hoist/ suspension to be used.
- Make sure sling does not have visible damage or frays.

Periodic inspection

Periodic inspection of the sling should be undertaken at least every 6 months. More frequent inspections may be required where the sling is used or cleaned more frequently than normal

Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings


Make sure not to damage or remove labels when cleaning.

Size and weight guide

Size	User weight - kg (lbs)
L	160-350 (353-772)
XL	300-500 (661-1102)

Sling checkpoints


1. The suspension is not too close to the users head
2. Legstraps are correctly attached according to the patients needs. Central legstraps crossed or legstraps on the same side.
3. The leg support is placed under users thigh without creases, and without risk of patient slipping out of the sling.
4. Place the sling in the level of the tailbone/ coccyx.

Sitting position - Sling application


1. Raise the back of the bed for a sitting position. Ask the user to lean forward, the patient may use the Immedia Bedstring to hold the position while sling is applied. Slide the sling down to the level of the tailbone/coccyx.

2. Pull the leg supports forward and position them under the user's thighs. Make sure the top of the sling is placed so it gives support to the users shoulders and neck.

3. Make sure the lifting straps are the same length on both sides and that the sling is placed under the users thighs without any creases.


4. Make sure the matching lifting straps are fastened to an equal length on left and right side and the sling supporting the body correctly without risk of patient sliding out of the sling. Make sure the user is in a comfortable and secure position in the sling.

Remove the sling in reverse order


Transfer

When moving the user, stand to the side of the person you are lifting. Make sure that arms and legs do not obstruct the seat, bed, etc.


1. Check that the sling is correctly fitted around the user and that the strap loops are correctly fitted to the suspension hooks.
2. Start lifting until the lifting straps are stretched without lifting the user. Ensure that all 6 loops of the sling are securely fastened to avoid the user slipping or falling out of the sling.
3. Lift user, and perform transfer.


Be careful during movement, the user may swing during turns, stops and starts. Be careful when manoeuvring close to furniture and similar to ensure the suspended user does not hit anything.

Not crossed

Crossing

Crossed

Depending on the patients needs, the central legstraps can be crossed or not. If crossed, pull one of the loops through the other. Connect all 6 lifting straps to the suspension.

Never hoist the user higher than necessary to carry out a lift.

Remove the sling gently to ensure the user does not roll over and fall off the bed.

Never leave a user unattended while being hoisted.


The hoist should not be used to lift or move users on sloping surfaces


Lying position - Sling application

We suggest the use of the Immedia SlingOn to apply the MediumBack PLUS in a correct and safe manner.


1. Fold the SlingOn in a way that makes it easy to slide and pull underneath the patient. The SlingOn must be double-layered with an opening at the top of the bed, this will make the application of the sling more efficient.

2. Slide the SlingOn under the patients head, a caregiver on each side of the patient should pull the SlingOn simultaneously down until it is fully unfolded and covering the patients body.



3. Slide the sling between the two layers of the SlingOn, starting from the top and working your way down.


4. Position the sling anatomically correct underneath the patient before removing the top layer of the SlingOn. Raise the back of the bed to make the patient have a sitting position before performing the lift.


5. Make sure the lifting straps are the same length on both sides and that the sling is placed under the users thighs without any creases. Depending on the patients needs, the central legstraps can be crossed or not. If crossed, pull one of the loops through the other. Connect all 6 lifting straps to the suspension.

Not crossed

Crossing

Crossed

Remove the sling in reverse order

Recycling

Recycling instructions are available at Etac.com

Troubleshooting

Symptom	Possible Cause/Action
The sling does not fit the user.	Wrong sling size. Try another size.
	The sling is too small when the leg supports do not reach halfway around the inner thigh and the suspension is too close to the user's face. Try another size.
	The sling is not correctly positioned. Pull the sling further down at the back when fitting it.
All the straps do not lift at the same time.	If the sling straps lift unevenly, it means that the sling can slide up at the back which can lead to the user sliding down into the opening of the sling. Make sure that the leg supports are positioned correctly on the thighs.
User is sliding through/out of the sling.	The sling is not connected correctly according to the users individual prescription.
	The user is wearing too slippery clothes. Change or remove some of the clothing.
	The sling is too big, try a smaller size or another sling

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